



Simply magnetic

FerroSelect™ Array



Closed, automated, and simple magnetic cell selection (coming soon)

We designed the FerroSelect™ Array to address many of the shortcomings of existing cell selection platforms. Developed from the ground up to take advantage of our proprietary array of neodymium iron boron rare earth magnets the instrument enables bulk collection of cells directly from either a quarter, or full leukopak.

The FerroSelect Array fully automated, closed, single-use cell selection platform brings high-volume cell capture to your benchtop. When used with our variety of positive and negative selection kits, cellular subsets are quickly and easily isolated from leukopaks.

SKU: 24-0003: FerroSelect Array

Research Use Only. Not for use in diagnostic procedures. Not for use in humans.

Key Features & Benefits



Single use consumable



Cell selection and depletion using validated protocols



User programmable with quick setup, operation, and cleanup




Figure 1: The FerroSelect Array consumable is designed to complement the capabilities of the platform, and let you quickly setup and execute a pre-programmed or customized selection procedure



Scan here

See how **FerroSelect** turns **simply magnetic** cell selection into **simply magic** results!

Simply scalable. Simply reproducible. Simply accessible. **Simply magnetic.**